


## Valnemulin

Code No.: **BIA-V1492**

Pack sizes: **5 mg, 25 mg**



Synonyms :

## Specifications

CAS #	: 101312-92-9
Molecular Formula	: C <sub>31</sub> H <sub>52</sub> N <sub>2</sub> O <sub>5</sub> S
Molecular Weight	: 564.8
Source	: <b>Semi-synthetic</b>
Appearance	: <b>White solid</b>
Purity	: <b>&gt;95% by HPLC</b>
Long Term Storage	: <b>-20°C</b>
Solubility	: <b>Soluble in ethanol, methanol, DMF or DMSO. Moderate water solubility.</b>

## Application Notes

Valnemulin is a semi-synthetic pleuromutilin prepared by sequential reaction of pleuromutilin tosylate with 2-amino-1,1-dimethylethylthiol, coupling the free amine with a protected valine and then deprotection to provide the final product. Valnemulin is a broad spectrum antibiotic used to control gastrointestinal infections in animals and shows no cross resistance to other antibiotic classes. Like all the pleuromutilins, valnemulin inhibits protein synthesis by binding to domain V of 23S RNA.

## References

1. Pleuromutilin derivatives a process for their preparation and their use. Berner H. et al., 1986, US Patent 4,675,330.