


Valnemulin hydrochloride

Code No.: **BIA-V1537**

Pack sizes: **25 mg, 100 mg**



Synonyms :

Specifications

CAS #	: 133868-46-9
Molecular Formula	: C ₃₁ H ₅₃ ClN ₂ O ₅ S
Molecular Weight	: 601.3
Source	: Semi-synthetic
Appearance	: White solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in ethanol, methanol, DMF or DMSO. Moderate water solubility.

Application Notes

Valnemulin hydrochloride is the salt prepared from valnemulin taking advantage of the basic amino group which protonates and readily forms the salt in hydrochloric acid solutions. The hydrochloride is the preferred formulation for pharmaceutical applications. Valnemulin is a broad spectrum antibiotic used to control gastrointestinal infections in animals, and shows no cross resistance to other antibiotic classes. Like all the pleuromutilins, valnemulin inhibits protein synthesis by binding to domain V of 23S RNA.

References

1. Pleuromutilin derivatives a process for their preparation and their use. Berner H. et al., 1986, US Patent 4,675,330.