

Small MWt unsaturated acids

Code No.: **BIA-AS6058**

Specifications

Sets are available as 1mg sets (1mg/vial) or 5mg sets (5mg/vial) of each of the following products:

Vial #	Compound	Code No.	CAS #	Mol. Formula	Mol. Wt.
1	Crotonic acid	BIA-C2438	3724-65-0	C ₄ H ₆ O ₂	86.1
2	Sorbic acid	BIA-S2437	110-44-1	C ₆ H ₈ O ₂	112.1
3	Tiglic acid	BIA-T2439	80-59-1	C ₅ H ₈ O ₂	100.1

Long Term Storage : **-20°C, protect from light**

Stability : **Stable for more than 1 year when stored at -20°C, protected from light**

Short Term Storage : **Stable at ambient temperature for 1-2 weeks, protected from light**

Shipping : **Ambient temperature**

Purity : **Minimum purity of >95% by HPLC**

Solubility :

Updated: 13 January 2023