

Valine homologues

Code No.: **BIA-AS6070**

Specifications

Sets are available as 1mg sets (1mg/vial) or 5mg sets (5mg/vial) of each of the following products:

Vial #	Compound	Code No.	CAS #	Mol. Formula	Mol. Wt.
1	D-Norvaline	BIA-N2168	2013-12-9	C ₅ H ₁₁ NO ₂	117.2
2	L-Norvaline	BIA-N2167	6600-40-4	C ₅ H ₁₁ NO ₂	117.2

Long Term Storage : **-20°C, protect from light**

Stability : **Stable for more than 1 year when stored at -20°C, protected from light**

Short Term Storage : **Stable at ambient temperature for 1-2 weeks, protected from light**

Shipping : **Ambient temperature**

Purity : **Minimum purity of >95% by HPLC**

Solubility :

Updated: 13 January 2023